

# Optogenetically Induced Spatiotemporal Gamma Oscillations in Visual Cortex

Fereshteh Arab<sup>1</sup>, Sareh Rostami<sup>2, 5</sup>, Mohammad Dehghani-Habibabadi<sup>3, 5</sup>, Vahid Salari<sup>1, 4 \*</sup>, and Mir-Shahram Safari<sup>2, \*</sup>

<sup>1</sup>Department of Physics, Isfahan University of Technology, Isfahan 84156-83111, Iran,

<sup>2, \*\*</sup> Neuroscience Research Center, Shahid Beheshti University of Medical Sciences, Tehran, 19615-1178, Iran,

<sup>3, \*\*</sup> Institute for Theoretical Physics, University of Bremen, 28359 Bremen, Germany,

<sup>4</sup>Department of Physical Chemistry, University of the Basque Country (EPV/EHU), Bilbao, Spain.

<sup>5, \*\*</sup> Equal contribution

<sup>1\*</sup> vahid.salari@ehu.eus

<sup>2\*</sup> safari@sbmu.ac.ir

## ABSTRACT

It has been hypothesized that Gamma cortical oscillations play important roles in numerous cognitive processes and may involve psychiatric conditions including anxiety, schizophrenia, and autism. Gamma rhythms are commonly observed in many brain regions during both waking and sleep states, yet their functions and mechanisms remain a matter of debate. Spatiotemporal Gamma oscillations can explain neuronal representation, computation, and the shaping of communication among cortical neurons, even neurological and neuropsychiatric disorders in neo-cortex. In this study, the neural network dynamics and spatiotemporal behavior in the cerebral cortex are examined during Gamma brain activity. We have directly observed the Gamma oscillations on visual processing as spatiotemporal waves induced by targeted optogenetics stimulation. We have experimentally demonstrated the constant optogenetics stimulation based on the ChR2 opsin under the control of the CaMKII $\alpha$  promotor, which can induce sustained narrowband Gamma oscillations in the visual cortex of rats during their comatose states. The injections of the viral vector [LentiVirus – CaMKII $\alpha$  – ChR2] was performed at two different depths, 200 and 500  $\mu$ m. Finally, we computationally analyze our results via Wilson-Cowan model.

## Introduction

One of the powerful experimental approaches in neuroscience to understand the spatiotemporal neuronal network is optogenetics. According to many types of research, optogenetics is a technique that involves using light to control cells in live tissues, typically neurons, that have been genetically modified to express light-sensitive ion channels. Optogenetics is a neuro-modulation method that uses a compound of techniques from optics and genetics to control and monitor individual neural activities in living tissues and precisely measure these manipulation effects in real-time,<sup>1-3</sup> and it enables us to manipulate neuronal activities in a millisecond. Optogenetics has furthered the fundamental scientific understanding of how specific cell types contribute to the function of biological tissues such as neuronal circuits in-vivo. Moreover, on the clinical side, optogenetic-driven research has led to insights into Parkinson's disease and other neurological and psychiatric disorders. Indeed, optogenetics has also provided different insights into neuronal codes relevant to Autism, Schizophrenia, drug abuse, anxiety, and depression.<sup>4-7</sup> Cortical oscillations are rhythmic patterns of neuronal activity that are synchronized across many neurons within or across brain regions that are classified based on frequency into Delta (0.5-4 Hz), Theta (4-8 Hz), Alpha (8-12 Hz), Beta (13-30 Hz), and Gamma (30-100 Hz) oscillations.<sup>8,9</sup> Cortical oscillation in Theta and Gamma frequency range has been hypothesized to play essential roles in numerous cognitive processes and may involve psychiatric conditions including anxiety, schizophrenia, and autism. Gamma rhythms are commonly observed in many brain regions during both waking and sleep states, yet their functions and mechanisms remain a matter of debate.

## History

Spatiotemporal Gamma oscillations can explain neuronal representation, computation, and the shaping of communication among cortical neurons, even neurological and neuropsychiatric disorders in neo-cortex.<sup>10-14</sup> Recent studies have used optogenetic stimulation to probe the mechanisms and neural circuits underlying Gamma oscillations' generation in the neo-cortex.<sup>15-19</sup> Moreover, different investigations have been done on optogenetically induced Gamma oscillations, such as compared Gamma oscillations in entorhinal-hippocampal system<sup>20</sup>, neuroprotective effects due to chronic gamma entrainment<sup>21</sup>, and a mechanistic

framework for modulation of pain by specific activity patterns in the primary somatosensory cortex<sup>22</sup>. In a recent optogenetic study, it has been elucidated the substrate oscillations in the prefrontal cortex via a protocol with a combination of optogenetics and electrophysiological recordings from neonatal mice, where the light was transferred to the layer II/III pyramidal neurons by in utero electroporation with Channelrhodopsin, boosting network oscillations within Gamma-Beta frequency range<sup>23</sup>. However, the activation of layer V/VI pyramidal neurons caused nonspecific network activation<sup>23</sup>. In an investigation<sup>24</sup>, it was expressed channelrhodopsin-2 under the CaMKII $\alpha$  promoter in mice and demonstrates that the cornuammonis area 1 (CA1) is capable of generating intrinsic Gamma oscillations in response to Theta input. This Gamma generator is independent of activity in the upstream regions, highlighting that CA1 can produce its own Gamma oscillation in addition to inheriting activity from the upstream regions, which supports the theory that Gamma oscillations predominantly function to achieve local synchrony and that a local Gamma generated in each area conducts the signal to the downstream region. Additionally, in another study, the optogenetic stimulation was performed on two macaque monkeys during their awake resting and reach and grasp states in the primary motor (M1) and ventral premotor (PMv) cortices of two subjects. They used the particular optogenetic construct including red-shifted opsin C1V1 under the control of the CaMKII $\alpha$  promoter and AAV5 viral vector.

### This study

In this paper, the effect of induced Gamma oscillations on visual power is investigated, and we check whether the visual activity of subjects increases these oscillations in the cortex or not. We would like to probe the connection between visual power and the power of induced Gamma oscillations in the second and the fifth layers of cortical. We investigate optogenetically induced Beta and Gamma oscillations in rat's visual cortex by Adeno Associated Virus (AAV) method. We express the channelrhodopsin2 in layers II and V, simultaneously with visual stimulation. As a novel aspect of our work, we have examined whether the visual response is affected by the optogenetic laser or not. In fact, the present study examines and analyzes the effect of induced Gamma oscillations by targeted optogenetics stimulation, recorded by intra-cortical one-channel extracellular techniques in rat's visual cortex during their comatose states. Injections of the viral vector [ LentiVirus - ChR2 – CaMKII $\alpha$ ] were performed at two different depths ( 200  $\mu$ m and 500  $\mu$ m). During trial stimulation, a 470-nm blue laser was delivered light through the polymer optical fiber. LFP trials in anesthetized rats in 4 states are recorded.

### Methodology

Three requirements are necessary to apply optogenetic techniques: selecting optogenetic molecular reagents (OMRs), targeted expression of OMRs in the neurons or regions of interest, and a light delivery system<sup>25</sup>. Neuronal control can be achieved using optogenetic actuators like Channelrhodopsin. An engineered viral vector is used to introduce the microbial opsin to a specific region of the organism, containing the optogenetics actuator gene attached to a recognizable promoter CaMKII $\alpha$ .<sup>3,26-30</sup> Various methods have been devised to express the OMRs in the cells to be controlled specifically, such as an adeno-associated virus (AAV) or lentivirus incorporating a target cell-specific promoter as a vector.<sup>5,31,32</sup> In this study, lentivirus was used. The changes in membrane potential, which show dynamic fluctuations in response to excitatory and inhibitory inputs, play key roles in excitable cells such as neurons and muscle cells. Electrophysiological approaches, including the micro-pipette method and the patch-clamp methods, allow recording fast voltage changes at a sub-millisecond level (e.g., action potentials), i.e. used here.<sup>33-35</sup>

From theoretical point of view, we model the cortex as a Wilson-Cowan neural field<sup>36</sup>. We consider laser stimulus as a source of an external current into the neuronal populations. During the laser onset, the outward current is applied to the neuronal population. We adjust the model parameters to explain our experimental results for describing variations in Beta and Gamma oscillations power.

### Computational Data Analysis

In this research, we use FIR1 filter with a band-pass type. First, we filter data in the range of 5 - 300 Hz, and then we focus on the range of Beta and Gamma frequencies. One of the quantities that obtain differences in four different conditions (Control, Visual, Laser, and Laser+Visual) is the mean power signal. Signal power in Gamma range from 30 - 150 Hz and Beta range 13 - 30 Hz were calculated for different trials. Equation 1 is the formula of discrete-time signal  $x(n)$  power.

$$\text{Signal Power} = \frac{1}{2N-1} \sum_{n=-N}^{n=N} |x(n)|^2 \quad (1)$$


To find the pure effect of Optogenetic laser at different depths, we define equation 2, as a biological response to laser stimulus.

$$\text{Biological Response} = (P_{\text{Laser+Visual}} - P_{\text{Laser}}) - (P_{\text{Visual}} - P_{\text{Control}}) \quad (2)$$

where P is the power of the signal.  $P_{\text{Visual}} - P_{\text{Control}}$  is any changes in the visual trials' power compared to the control trials. If the quantity gets positive, it means that the visual stimulus has a positive effect in Beta band frequencies.

$P_{\text{Visual+Laser}} - P_{\text{Laser}}$  shows changes caused by vision. If Biological Response is equal to zero, there is no effect from optogenetic stimulus on visual responses.

In this simulation, whenever the Laser is on, the corresponding LFPs data are not considered. We concentrate on the biological response between two pulses. As figure 1 shows, blue lines are non-zero at the time of Laser is switched on, and the corresponding LFP is shown in red color at the same time. Figure 1, shows at the time of Laser is switched on, there is a jump in LFP data, which gives a high frequency in the analysis. By removing these data and focusing on time between two pulses, artifacts do not produce high frequencies.



**Figure 1.** (Color online) Blue line shows amplitude of laser in time. Red line is corresponding LFP during laser stimulus.

### Modeling Approach

We use the Wilson Cowan model (as a standard model in optogenetics)<sup>37,38</sup> for modeling our study as follows in the equation 3.

$$\begin{aligned} \tau_e \frac{dE(r,t)}{dt} &= -E(r,t) + f(W_{ee} * E(r,t) + W_{ei} * I(r,t) + \theta_e + J_e) \\ \tau_i \frac{dI(r,t)}{dt} &= -I(r,t) + f(W_{ie} * E(r,t) + W_{ii} * I(r,t) + \theta_i + J_i) \end{aligned} \quad (3)$$

$E$  and  $I$  represent mean firing rates of the excitatory and inhibitory population at the point  $r$ ,  $(x,y)$ , respectively. Each population has a specific membrane time constant  $\tau$ , and  $\theta$  is a continuously synaptic current, visual stimulation, or other currents to the cells. The variables  $J_e$  and  $J_i$  are excitatory and inhibitory currents to each cell respectively, due to light in the optogenetics setup<sup>39</sup>. The Gaussian distribution is considered for the domain and efficacy of inhibitory and excitatory synapses, and the parameter  $W$  is defined based on this distribution, i.e. the synaptic influence of each population on itself and other cells, as follows


$$W = \frac{w}{\sqrt{2\pi}s} \exp\left(-\frac{x^2}{2s^2}\right) \quad (4)$$

$w$  and  $s$  depend on inhibitory and excitatory populations are different. Figure 2 b shows inhibitory and excitatory cells power and effect's domain. There are four types of connections, and  $w$  is changed depending on the type of connection. inhibitory to excitatory,  $w_{ie}$ , inhibitory to inhibitory,  $w_{ii}$ , excitatory to excitatory,  $w_{ee}$ , excitatory to inhibitory,  $w_{ei}$ . In our study  $w_{ii}$  and  $w_{ie}$  is  $-0.4$  and  $w_{ee}$  and  $w_{ei}$  is  $+0.1$ .  $s_i$  is  $0.04$  and  $s_e$  is  $0.02$ .

$$K_{ie}(x) * I(x,t) = \int K_{ie}(x-y) I(x,t) dy \quad (5)$$

In this study, we have considered 100 inhibitory and 100 excitatory cells, and each cell dose not have connection to itself. The function  $f$  is the threshold function, which is written in equation 5.

$$f(x) = \frac{1}{1 + e^{-4(u-1)}} \quad (6)$$



**Figure 2.** (color online) Left: Two different types of inputs to the population. Right: Excitatory and inhibitory population domain and efficacy.

where  $u$  is local synaptic input at spatial position  $r$  at time  $t$ . LFP is a contribution of excitatory and inhibitory populations activity, equation 6.

$$\text{LFP} = 0.8E + 0.2I \quad (7)$$

Generating oscillations depend on the inputs to the populations, and there are different inputs to the populations in different trials. We consider that there is no input to inhibitory cells, therefore,  $\theta_i = 0$  and  $J_i = 0$ . Figure 2 a is shown different type of current to excitatory population. A constant input plus a periodic Beta frequency, 15 Hz, input  $1 + \sin(15\pi t)$  to the excitatory population ensures Beta and Gamma oscillations in the control state. in the beginning  $\tau_i = 4ms$  and  $\tau_e = 2ms$  We consider a periodic triangular laser stimulus with an absolute amplitude to the excitatory neurons, and laser time onset is the same as the experiment. Opsin's main effect is on excitatory neurons. Therefore laser stimulus is considered only for the excitatory cells.

## Materials and Methods

The general setup of the experiment is shown in Figure 3. The details of experiment are as follows:

### Experimental Task

#### Animal preparations and virus injection


Ten naive rats of either sex, weight 120 – 150 g, were kept under standard conditions at  $20 \pm 2$  °C and light: dark hours of 12:12 with free access to water and standard food. All experimental procedures were performed according to the Guide for the Care and Use of Laboratory Animals (National Institute of Health Publication No. 80-23, revised 1996) and were approved by the ethics committee, Shahid Beheshti University of Medical Sciences, ethical code IR.SBMU.PHNS.REC. 1398.024. Animals were anesthetized by ketamine/xylazine (ip, 80/2.5 mg/kg). Their body temperature maintained at  $36.0 - 36.5$ °C by a rectal thermos probe feeding back to a heating pad (ATC-402, Unique Medical). A small custom-made head plate to stabilize the animal's head under virus injection and LFP recording attached over the occipital region of the left hemisphere. After determining of cranial window over primary visual cortex (V1) (3-9 mm posterior to Bregma; 1-5 mm lateral), a small part of the skull and dura mater were removed and exposed cortex was covered by Ringer's solution. We injected viral vector encoding Channelrhodopsin-2 under the control of CaMKII promotor (1.5- 2  $\mu$ l) in layers II-IV (200 - 500  $\mu$ m) of V1 by applying positive pressure to glass electrode (pipette) preparing for LFP recording in same area two weeks later.

#### Data acquisition and visual stimulation

Square-wave gratings (0.08 cycle/deg) at  $100 \text{ mW/mm}^2$  contrast was moved on an LCD monitor in 12 directions at 30 degree steps (0-360 degrees). These 12 patterns of visual stimuli (2 second before, 3 second during and 2 second after the presentation of each stimulus, respectively) were presented 3 times in a randomly shuffled order. The LCD monitor covered  $80 \times 50$  degrees of the visual field at a viewing distance of 28 cm.

#### Optogenetic stimulation

A blue laser light (470 nm) was directed into an optical fiber (700  $\mu$ m diameter), mounted on a holder for positioning in the recording chamber. The laser's output function was tested with the optical power meter at the tip of the optical fiber. Only a low total output power may be needed to achieve ChR2 activation. We used three trials in different laser intensities (input:  $50 \text{ mW/mm}^2$  and output:  $5 \text{ mW/mm}^2$ ), (input:  $75 \text{ mW/mm}^2$  and output:  $31 \text{ mW/mm}^2$ ), (input:  $100 \text{ mW/mm}^2$  and output  $80 \text{ mW/mm}^2$ ), respectively. Each trial consists 10 pulses at 20 Hz (2 s light off, 3 s light on (5 ms delay, 45 ms during)).



**Figure 3.** Setup of the experiment. Here, in order to record V1 data in the left hemisphere, the monitor was located in the right side of animal instead of front side.

### LFP recording

Animal were anesthetized with Urethane (ip, 1.5 mg/g body weight) and then placed in a stereotaxic frame (craniotomy). LFP recording from layer II (200  $\mu\text{m}$ ) and layer IV (500  $\mu\text{m}$ ) of V1 were performed by recording electrodes pulled from borosilicate glass capillary with filaments (0.86 mm inner diameter, 1.5 mm outer diameter) filled by Ringer's solution. The resistance of these electrodes was 500 k $\Omega$ . LFPs were recorded as four groups (three trials for each): 1) without visual stimulation, 2) without visual stimulation but with optogenetic stimulation, 3) with visual stimulation, 4) with visual stimulation and optogenetic stimulation. Using a molecular device amplifier (Axopatch 200B, Molecular Devices), recorded data was sampled at 20 kHz, filtered at 2-5 KHz, digitized at 20 KHz, and fed into a personal computer with an NI-DAQ board (PCI-MIO-16E-4, National Instruments).


### Results

Figure 4 shows Local Field Potential (LFP) sample in Laser+Visual stimulus. In the two first seconds there is no laser stimulus, and it is switched-on after that. We first focus on variations in Gamma and Beta frequencies power in different conditions.


Figure 5a shows Gamma power in different conditions. As shown in figure 5, Gamma power increases in laser trials in both 200 and 500  $\mu\text{m}$ , which means optogenetic stimulus induces more powerful Gamma oscillations in comparison with visual trials. However, optogenetic laser leads to Gamma oscillations; in 200  $\mu\text{m}$ , Gamma power is more significant than 500  $\mu\text{m}$ . When the subject is under both visual and laser stimuli, Gamma power decreases in both depths (except 500  $\mu\text{m}$  with 100  $\text{mW}/\text{mm}^2$  intensity laser); it starts to increase when lasers' intensity increases.

Figure 5b shows Beta power in different conditions. When there is only laser stimulus, Beta power decreases in both depths. When the subject is under both visual and laser stimuli, Beta power decreases in 200  $\mu\text{m}$  (except in 100  $\text{mW}/\text{mm}^2$  intensity laser). In 500  $\mu\text{m}$  in Laser+Visual trials, the first Beta power increases and starts to decrease by implementing high laser intensity. In order to see if optogenetic leads to have different conditions in vision we compare  $P_{\text{Visual}} - P_{\text{Control}}$  and  $P_{\text{Laser+Visual}} - P_{\text{Laser}}$ , we define it as biological response (see materials and methods). The results for 500  $\mu\text{m}$  and 200  $\mu\text{m}$  are summarized in table 1 and table 2 respectively. Both tables show that the Optogenetic laser has a different influence on visual power in both depths. The laser's impact is mainly negative in 200  $\mu\text{m}$  for both Beta and Gamma oscillations, while for laser intensities 50  $\text{mW}/\text{mm}^2$  and 100  $\text{mW}/\text{mm}^2$  are positive for both Beta and Gamma oscillations at 500  $\mu\text{m}$ .

Our approach to understanding the basic behind the results is Wilson-Cowan model<sup>36</sup>. The specific range of injected currents leads to Gamma oscillations in population activity<sup>40</sup>. Besides, Optogenetic stimulus active inhibitory and excitatory neurons by changing their firing rates<sup>36</sup>. Because of opsin's type, the Optogenetic laser excites excitatory neurons; therefore,



**Figure 4.** LFP in Laser+Visual condition in different depths.



**Figure 5.** (Color online) a: Mean Gamma power; b: Mean Beta power in Laser+Visual condition and different depths. Red for 500 μm and blue for 200 μm.

	$P_{\text{Visual}} - P_{\text{Control}}$	$P_{\text{Visual+Laser}} - P_{\text{Laser}}$			$(P_{\text{Visual+Laser}} - P_{\text{Laser}}) - (P_{\text{Visual}} - P_{\text{Control}})$		
laser intensity mW/mm <sup>2</sup> :	0	50	75	100	50	75	100
Frequency range: Beta	0.02±0.0087	0.08±0.0188	0.02±0.0155	-0.04±0.0187	0.06±0.0275	±0.0242	-0.06±0.0247
Gamma	±0.0069	-0.02±0.0147	-0.01±0.0169	0.14±0.0114	-0.02±0.0216	-0.01±0.0238	+0.14±0.0183

**Table 1.** Biological response in 500  $\mu\text{m}$ .

	$P_{\text{Visual}} - P_{\text{Control}}$	$P_{\text{Visual+Laser}} - P_{\text{Laser}}$			$(P_{\text{Visual+Laser}} - P_{\text{Laser}}) - (P_{\text{Visual}} - P_{\text{Control}})$		
laser intensity mW/mm <sup>2</sup> :	0	50	75	100	50	75	100
Frequency range: Beta	0.02±0.0113	-0.06±0.0128	-0.13±0.0131	-0.04±0.024	-0.08±0.0241	-0.15±0.0244	-0.06±0.0353
Gamma	0.02±0.0089	-0.25±0.0196	-0.02±0.0208	0.02±0.0165	-0.027±0.0285	-0.04±0.0297	±0.0254

**Table 2.** Biological response in 200  $\mu\text{m}$ .

we consider a pulse input to the excitatory neurons; however, this current does not significantly change the Gamma power. Figure 6 shows the amount of Gamma power per power of the whole signal versus the different excitatory time constants. Figure 6 shows while excitatory time constant is decreasing, the power of Gamma is increasing. It can describe why the Gamma power in 200 and 500  $\mu\text{m}$  increases in laser mode compared to visual and control modes. Therefore we conclude the Gamma power changes might arise due to the changes in population time constants.

We concluded in Laser+Visual states that vision acts as a resistance and prevents the excitatory populations from decreasing time constant. As shown in figure 5, Gamma power in Laser+Visual is less than in laser trials; however, it is still more than visual and control states.

As figure 6 b shows, Gamma oscillations show a monotonic behavior; however, Beta oscillations are nonmonotonic (Figure 6 ). Figure 6 shows the nonmonotonic behavior of Beta oscillations in both depths.


Therefore, our modeling approach suggests that Beta oscillations can increase or decrease depending on initial conditions; figure 6 shows this nonmonotonic behavior in empirical results. As described before, Optogenetics is amplifying Gamma oscillations by decreasing the time constant of the excitatory population. If the initial time constant of excitatory cells was near the peak in figure 6, then the decreasing time constant leads to less Beta oscillation, which is in agreement with Figure 6 in laser trials. However, in the laser+visual trials, while vision acts as a resistance in 500  $\mu\text{m}$ , it is a resonator in 200  $\mu\text{m}$  and decreases Beta oscillations significantly. One can say, vision tries to push time constant to peak in 500  $\mu\text{m}$  while away from the maximum in 200  $\mu\text{m}$ .

## Discussion

Gamma and Beta Oscillation's average power in control and visual trials at 500  $\mu\text{m}$  is slightly more than 200  $\mu\text{m}$ . However, Gamma oscillations are amplified at 200  $\mu\text{m}$ , and Beta oscillations (except for laser trials at the intensities of 50 and 75 mW/mm<sup>2</sup>) show similar behavior at 500  $\mu\text{m}$ . Gamma oscillations at all trials at 200  $\mu\text{m}$  have more power than 500  $\mu\text{m}$ . Except for laser trials at intensities of 50 and 75 mW/mm<sup>2</sup>, Beta oscillations power at 500  $\mu\text{m}$  is slightly more than 200  $\mu\text{m}$ . Our modeling approach (see section "Modeling Approach") suggests that results at 50 and 75 mW/mm<sup>2</sup> in laser trials arise from Beta oscillations' nonmonotonic behavior versus the changing time constant of excitatory populations.

### Layer II

The Gamma power in all three laser intensities is higher than in the control mode, where our model describes this behavior by shifting the excitatory time constant to the higher amounts. Moreover, Gamma power in all three Laser+Visual trials is higher than the visual trials. Therefore, our results indicate a positive effect of the Optogenetics on the induction of Gamma oscillations. The power of Gamma oscillation in laser trials is higher than in Laser+Visual trials. While our modeling approach suggests Optogenetics laser shift excitatory time constant to the higher amounts; it seems vision acts as resistance and prevents shifting a



**Figure 6.** Gamma and Beta power versus  $\tau_e$ : monotonic behavior of Gamma power, and nonmonotonic behavior of Beta power by considering  $\tau_i = 0.004$ .

lot. The Beta oscillations power in laser and Laser+Visual trials are less than control and visual trials, and Optogenetics laser is reducing the power in all trials. However, the Optogenetics laser reduces Beta power; Laser+Visual trials have less Beta power than laser trials. Our modeling approach suggests that it is due to Beta oscillations' nonmonotonic behavior, and we predict vision and Optogenetics laser cooperate to reduce Beta oscillations power in  $200\mu\text{m}$ . In contrast, as discussed above, vision makes it resistant to prevent sharp increases in Gamma oscillations. A comparison of two laser and Laser+Visual trials indicates that the laser trials have more substantial beta oscillation at  $50$  and  $75\text{ mW/mm}^2$  intensities. However, at the  $100\text{ mW/mm}^2$  intensity, the Laser+Visual state's beta-oscillation is more than laser trials. The excessive intensity of the laser seems to reverse the behavior entirely.

### Layer V

The laser trials' Gamma power is higher than the control trials, and the Laser+Visual trials are higher than the visual trials. Therefore Optogenetic amplifies the power of Gamma oscillations at  $500\mu\text{m}$ . By considering laser and Laser+Visual trials, vision prevents sharp increment in Gamma power; even with different laser intensities. Beta oscillations power in laser trials at  $50$  and  $75\text{ mW/mm}^2$  compared to the control trials is reducing. However, both laser and control are approximately at the same level of Beta power in intensity. Our modeling approach (Figure 4) suggests that if the initial time constant of excitatory is near the peak, changing the time constant can decrease Beta power. Increasing Gamma power while decreasing Beta power is due to Gamma oscillations' monotonic behavior and nonmonotonic behavior of Beta oscillation.

Here we observed two different vision behavior in  $200\mu\text{m}$  and  $500\mu\text{m}$ . While vision in cooperation of laser reduced the Beta power in  $200\mu\text{m}$ , it shows the opposite behavior in  $500\mu\text{m}$ . According to our modeling approach, it could be due to the different time constant of populations in different depths. Generally, the Optogenetic laser amplifies Beta oscillations' power and reduces the Beta oscillations' power at both depths. Biological responses at  $200\mu\text{m}$  are negative in all cases, and at  $500\mu\text{m}$ , it is highly dependent on laser power. Laser intensity and visual power have an indirect relationship in Beta oscillations; while it is a direct relationship for Gamma oscillations, we have concluded cortical responses for Gamma and Beta oscillations are opposite.

In our experiment, the opsin affects the excitatory population, and then we consider a laser pulse current to the excitatory population. The different laser injection amplitudes can make the excitatory population faster or slower via changing the time constant. As the excitatory population gets fast, then the inhibitory population will get fast too, and therefore some oscillatory behaviors could emerge.


We conclude that when optogenetic stimulus changes the time constant of the inhibitory and excitatory neurons, specific amounts of injected current increase or decrease the frequencies power.

### Conclusion

Our simulations propose that considering different time constant of each population plays the most crucial role in modeling the stimulus<sup>36</sup>. The proposed model suggests that each population's time constant is changed under the optogenetic stimulus's influence, and the ratio of the time constants of the inhibitory and excitatory population explains the Optogenetic behavior. The laser excites the excitatory cells, and it significantly changes Gamma power by changing the time constant of excitatory. Figure 7 a shows no significant change in Gamma power in the whole time. There is no laser stimulus in the one and a half seconds,



while there is a laser stimulus in the rest. Figure 7a agrees with Laser+Visual trial in 200 microns, Figure 7c. Figure 7b shows a significant change in Gamma power when the time constant of excitatory neurons are more than inhibitory neurons. Figure 7d can describe the data in figure 7d which is a trial in 200 microns in Laser+Visual trial. As a conclusion, optogenetic



**Figure 7.** (Color online) Variation of Gamma power in response to different amounts of excitatory and inhibitory time constants. a:  $\tau_e < \tau_i$  b:  $\tau_e > \tau_i$  c: The fifth trial in Laser Visual in 200 microns d: The Third trial in Laser+Visual in 500 microns.

stimulation in both depths of visual cortex induces gamma oscillations for all applied laser intensities, while beta oscillations are diminished for some intensities. In fact, the responses to Beta and Gamma oscillations are different; it is monotonic in Gamma oscillations but non-monotonic in Beta oscillations.

## Acknowledgment

All authors gratefully acknowledge financial support from their affiliated institution and thank Mehdi Aslani for graphical representation of the experimental setup. MDH is very thankful for valuable discussions with Udo Ernst for the theoretical modelling.

## References

1. Deisseroth, K. *et al.* Next-generation optical technologies for illuminating genetically targeted brain circuits (2006).
2. March, J. A. History of optogenetics .
3. Yizhar, O., Fenno, L. E., Davidson, T. J., Mogri, M. & Deisseroth, K. Optogenetics in neural systems. *Neuron* **71**, 9–34 (2011).
4. Deisseroth, K. Optogenetics. *Nature Methods* **8**, 26 (2011).
5. Fenno, L., Yizhar, O. & Deisseroth, K. The development and application of optogenetics. *Annual Review of Neuroscience* **34** (2011).
6. Yizhar, O., Fenno, L. E., Davidson, T. J., Mogri, M. & Deisseroth, K. Optogenetics in neural systems. *Neuron* **71**, 9–34 (2011).
7. Deisseroth, K. Optogenetics and psychiatry: applications, challenges, and opportunities. *Biological Psychiatry* **71**, 1030–1032 (2012).

8. Hari, R. & Salmelin, R. Human cortical oscillations: a neuromagnetic view through the skull. *Trends in Neurosciences* **20**, 44–49 (1997).
9. Buzsáki, G. & Draguhn, A. Neuronal oscillations in cortical networks. *Science* **304**, 1926–1929 (2004).
10. Buzsáki, G. & Wang, X.-J. Mechanisms of gamma oscillations. *Annual Review of Neuroscience* **35**, 203–225 (2012).
11. Tiesinga, P. & Sejnowski, T. J. Cortical enlightenment: are attentional gamma oscillations driven by ing or ping? *Neuron* **63**, 727–732 (2009).
12. Ray, S. & Maunsell, J. H. Do gamma oscillations play a role in cerebral cortex? *Trends in Cognitive Sciences* **19**, 78–85 (2015).
13. Womelsdorf, T., Fries, P., Mitra, P. P. & Desimone, R. Gamma-band synchronization in visual cortex predicts speed of change detection. *Nature* **439**, 733 (2006).
14. Van Kerkoerle, T. *et al.* Alpha and gamma oscillations characterize feedback and feedforward processing in monkey visual cortex. *Proceedings of the National Academy of Sciences* **111**, 14332–14341 (2014).
15. Sohal, V. S. Insights into cortical oscillations arising from optogenetic studies. *Biological Psychiatry* **71**, 1039–1045 (2012).
16. Csicsvari, J., Jamieson, B., Wise, K. D. & Buzsáki, G. Mechanisms of gamma oscillations in the hippocampus of the behaving rat. *Neuron* **37**, 311–322 (2003).
17. Traub, R., Whittington, M., Colling, S., Buzsáki, G. & Jefferys, J. Analysis of gamma rhythms in the rat hippocampus in vitro and in vivo. *The Journal of Physiology* **493**, 471–484 (1996).
18. Yamamoto, J., Suh, J., Takeuchi, D. & Tonegawa, S. Successful execution of working memory linked to synchronized high-frequency gamma oscillations. *Cell* **157**, 845–857 (2014).
19. Buhl, E. H., Tamás, G. & Fisahn, A. Cholinergic activation and tonic excitation induce persistent gamma oscillations in mouse somatosensory cortex in vitro. *The Journal of Physiology* **513**, 117–126 (1998).
20. Butler, J. L., Hay, Y. A. & Paulsen, O. Comparison of three gamma oscillations in the mouse entorhinal–hippocampal system. *European Journal of Neuroscience* **48**, 2795–2806 (2018).
21. Adaikkan, C. & Tsai, L.-H. Gamma entrainment: Impact on neurocircuits, glia, and therapeutic opportunities. *Trends in neurosciences* **43**, 24–41 (2020).
22. Tan, L. L. *et al.* Gamma oscillations in somatosensory cortex recruit prefrontal and descending serotonergic pathways in aversion and nociception. *Nature communications* **10**, 1–17 (2019).
23. Bitzenhofer, S. H. *et al.* Layer-specific optogenetic activation of pyramidal neurons causes beta–gamma entrainment of neonatal networks. *Nature communications* **8**, 1–13 (2017).
24. Butler, J. L., Mendonça, P. R., Robinson, H. P. & Paulsen, O. Intrinsic cornu ammonis area 1 theta-nested gamma oscillations induced by optogenetic theta frequency stimulation. *Journal of Neuroscience* **36**, 4155–4169 (2016).
25. Yawo, H., Asano, T., Sakai, S. & Ishizuka, T. Optogenetic manipulation of neural and non-neural functions. *Development, growth & differentiation* **55**, 474–490 (2013).
26. Zhang, F. *et al.* Red-shifted optogenetic excitation: a tool for fast neural control derived from *volvox carteri*. *Nature Neuroscience* **11**, 631 (2008).
27. Nagel, G. *et al.* Channelrhodopsin-1: a light-gated proton channel in green algae. *Science* **296**, 2395–2398 (2002).
28. Berthold, P. *et al.* Channelrhodopsin-1 initiates phototaxis and photophobic responses in *chlamydomonas* by immediate light-induced depolarization. *The Plant Cell* **20**, 1665–1677 (2008).
29. Mattis, J. *et al.* Principles for applying optogenetic tools derived from direct comparative analysis of microbial opsins. *Nature Methods* **9**, 159 (2012).
30. Zhang, F. *et al.* The microbial opsin family of optogenetic tools. *Cell* **147**, 1446–1457 (2011).
31. Zhang, F., Wang, L.-P., Boyden, E. S. & Deisseroth, K. Channelrhodopsin-2 and optical control of excitable cells. *Nature Methods* **3**, 785 (2006).
32. Gradinaru, V. *et al.* Molecular and cellular approaches for diversifying and extending optogenetics. *Cell* **141**, 154–165 (2010).
33. Scanziani, M. & Häusser, M. Electrophysiology in the age of light. *Nature* **461**, 930 (2009).

34. Ting, J. T., Daigle, T. L., Chen, Q. & Feng, G. Acute brain slice methods for adult and aging animals: application of targeted patch clamp analysis and optogenetics. In *Patch-Clamp Methods and Protocols*, 221–242 (Springer, 2014).
35. Wang, G. *et al.* An optogenetics-and imaging-assisted simultaneous multiple patch-clamp recording system for decoding complex neural circuits. *Nature Protocols* **10**, 397 (2015).
36. Heitmann, S., Rule, M., Truccolo, W. & Ermentrout, B. Optogenetic stimulation shifts the excitability of cerebral cortex from type i to type ii: oscillation onset and wave propagation. *PLoS Computational Biology* **13**, e1005349 (2017).
37. Wilson, H. R. & Cowan, J. D. A mathematical theory of the functional dynamics of cortical and thalamic nervous tissue. *Kybernetik* **13**, 55–80 (1973).
38. Wilson, H. R. & Cowan, J. D. Excitatory and inhibitory interactions in localized populations of model neurons. *Biophysical Journal* **12**, 1–24 (1972).
39. Lu, Y. *et al.* Optogenetically induced spatiotemporal gamma oscillations and neuronal spiking activity in primate motor cortex. *Journal of Neurophysiology* **113**, 3574–3587 (2015).
40. Onslow, A. C., Jones, M. W. & Bogacz, R. A canonical circuit for generating phase-amplitude coupling. *PLoS One* **9**, e102591 (2014).